



User manual




Technical Specifications

Model	AL300
LED Type	5W LED
LED QTY	60pcs
LED power	300w
Wavelength	630nm/660nm/810nm/830nm/850nm(Standard)
Beam Angle	30 degree
Cooling Fans Quantity	1 pc
Product size	316*214*60mm
EMF Emission	0 μ T @ 3inches
Sound(dB)@15cm	52 DB
Flicker	No Flicker
Input voltage	AC 110-240V 50/60Hz
Working Temperature	- 20°C to 40°C
Life Span	100000 hours
Warranty	2 Years

Technical Specifications

Model	AL300
LED Type	5W LED
LED QTY	60pcs
LED power	300w
Wavelength	630nm/660nm/610nm/830nm/850nm(Standard)
Beam Angle	30 degree
Cooling Fans Quantity	1 pc
Product size	316*214*60mm
EMF Emission	0uT @ 3inches
Sound(dB)@15cm	52 DB
Flicker	No Flicker
Input voltage	AC 110-240V 50/60Hz
Working Temperature	- 20°C to 40°C
Life Span	100000 hours
Warranty	2 Years




IR



RED



Power Switch:






Near Infrared light at 810 nm
830 nm 850 nm

Near Infrared light penetrates deep
Your body-enhancing muscle recovery
reducing joint pain

Product features:

This type of light will go the deepest into your skin. It will help to address the clarity and tone of your skin as well as the deeper bone and muscle. Great for elastin and collagen production ,



Red Light at 630 nm 660 nm

Most of red light absorbed by your skin, leading to increase collagen production and optimal skin health

Near Infrared light at 810 nm 830 nm 850 nm

Near Infrared light penetrates deeper within your body- enhancing muscle recovery and reducing joint pain

Infrared 750-1000 nm

Results Boosted: Penetrates deeper into the skin for amplified results.

Red

625-700 nm

Anti-Aging: Combats wrinkles and increases collagen production for smoother, younger-looking skin.

Medical Indications (many)

- Pain relief
- Increases
- Sprains Enzymes boosted to improve
- For healing and energy Boosts DNA
- Healing Alleviates allergy symptoms
- Speeds healing of burns
- Skin repair

- Boosts liver function
- Beneficial for lungs
- Promotes circulation
- Beneficial for bones
- Beneficial for heart
- Beneficial for kidney
- Promotes wound he

HOW TO USE IT.

1. Distance: Apply the light from 3-60 cm to the body

2. Session Time: Exposure each area for 5-30 minutes, depen

Repeat on other areas if desired

3. Frequency: Do this 2 times to 14 times a week until results a

4. Maintenance: Continue using several times a week

Near Infrared (NIR)

- Pain relief
- Increases
 - Speeds wound healing
 - Healing of muscle and ligament
- Sprains Enzymes boosted to improve body
 - Processes Boosts mitochondria function in cells
- For healing and energy Boosts DNA in body cells for better
- Healing Alleviates allergy symptoms
- Speeds healing of burns
- Skin repair

Red

- Boosts liver function
- Beneficial for lungs
 - Skin repair
 - Reduces scar tissue
- Promotes circulation Pain relief
 - Enhances the senses
 - Improves nutritional metabolism
- Beneficial for bones
 - Regenerates blood
- Beneficial for heart
- Beneficial for kidneys Boosts liver function
 - Promotes muscle healing
- Promotes wound healing

How to use it:

1. Distance: Apply the light from 3-60 cm to the body
2. Session Time: Exposure each area for 5-30 minutes, depending on distance.
Repeat on other areas if desired
3. Frequency: Do this 2 times to 14 times a week until results are felt
4. Maintenance: Continue using several times a week

Tips:

This light is for indoor use only, and please do not use in the wet environment. Please cut off the power when stop working. Please do not place it near any heat source.

Can this hurt my eyes?

The typical human eye responds to wavelengths 400nm to 700nm. The light coming from this device is either visible red or invisible infrared.

Infrared light is beneficial to the eyes, just as it is elsewhere on the body. Safety goggles are not necessary to use with infrared.

Red light is bright and can be uncomfortable to use directly on the eyes. Closing your eyes is sufficient to use red light directly on the face.

Tips:

This light is for indoor use only, and please do not use in the wet environment.
Please cut off the power when stop working.
Please do not place it near any heat source.

Can this hurt my eyes?

The typical human eye responds to wavelengths 400nm to 700nm. The light coming from this device is either visible red or invisible infrared.
Infrared light is beneficial to the eyes, just as it is elsewhere on the body. Safety goggles are not necessary to use with infrared.
Red light is bright and can be uncomfortable to use directly on the eyes. Closing your eyes is usually enough to use red light directly on the face.

CAUTION AND WARNING :

1. The product is not allow to expose to high temperature and high humidity for long time.
2. Lamps must be used at predetermined voltage and frequency.
3. The metal parts of lamps should not be polished with powder polisher, the dust behind the lamps should be cleaned with dry cloth or duster, which can be cleaned once a week.