

## Certification

Execution standard: GB4706.8-2008

Product name: Infrared Blanket

Test: 

Product date: 20



## Warranty card

### Records of the warranty

| Warranty date | Failure and processing methods | Completion Date | Customer signature |
|---------------|--------------------------------|-----------------|--------------------|
|               |                                |                 |                    |
|               |                                |                 |                    |
|               |                                |                 |                    |

### user information

|                  |  |
|------------------|--|
| Product model    |  |
| Product number   |  |
| Date of purchase |  |
| Customer name    |  |

## INSTRUCTIVE MANUAL

### FUERLE INFRARED BLANKET

**IMPORTANT:** Please read these instructions before using your FUERLE infrared blanket, Keep this manual handy.

#### Manufacturer :

Guangzhou Fuerle Health Technology Co., Ltd

Building No.1 Phase No.3 , Zhongheng International Commercial Centre

Yingbing avenue, Huadu district

Guangzhou City,

Guangdong,China

CB UK CE RoHS 

**MADE IN CHINA**

**Nothing is more important than health**

...



## READ BEFORE USE

### IMPORTANT:

1. Before use, please measure the blood pressure to be normal, remove the metal objects from your body, and it is more conducive to penetration if you enter the cabin naked (or wear disposable underwear), with your head exposed.

2. Before, during and after use, please drink 200-300 ml warm boiled water (can add a little salt) every time, to replenish water in the body, which is conducive to flush toxins out of the body.

3. For the first time, please use the low or medium temperature range, and the recommended time is 30-45 minutes.

4. If conditions permit, please use it after bathing.

5. Don't use it one hour before and after meals. Don't use it if you are too hungry or full. Don't use it if you are extremely weak.

6. Do not blow the air conditioner after use, and do not take a bath within 2 hours after use.

7. It is forbidden to be used by people with high blood pressure. It is forbidden to be used by people whose heat is not sensitive, whose heart is equipped with pacemaker, and who can't take care of themselves.

8. It is forbidden to use it for pregnant women, infants and children.

## PRECAUTIONS:

1. Must be used in the area where the leakage protection switch is installed.
2. must be connected to a socket with a neutral wire (ground wire) 6A.
3. It is strictly forbidden to use it during thunderstorms and lightning.
4. High temperature state (above 60°C), can not continuously exceed 90 minutes.
5. It is strictly forbidden to use it when the appearance of the product is damaged, the power plug is damaged, or the temperature controller is faulty.
6. It is strictly forbidden to use the product when folded.
7. After use, only the front side can be folded, and the back side can not be folded. (The side with the photon lamp can be folded inside)
8. It is strictly forbidden to use it together with cosmetics or drugs with strong acid for weight loss.
9. Don't use it naked or wearing thick coats.
10. Pull out the plug and cut off the power when leaving or not in use.



1. While in use, do not place anything (including the controller) on any hot surface of the blanket.
2. Do not fold the blanket while it is hot and cooling down. This can cause damage to the heating wires and the interior material of the blanket. Let it cool down completely before folding.

## Operating instruction

1. Connect the controller and connect the power.
2. Press the power button.
3. Press the start button.
4. Press temperature ▲ adjust temperature.
5. Press time ▲ adjust the use time.
6. After use, please unplug the power.



## Operating instruction

1. Connect the controller and connect the power.
2. Press the power button.
3. Press the photon button to turn on the photon function.
4. Press temperature ▲ ▼ adjust temperature.
5. Press time ▲ ▼ adjust the use time.
6. After use, please unplug the power.



Thank you for your choice Fuerle.

Fuerle infrared sauna blankets provide a long list of potential benefits-from detoxifying the body,Buring calories,to rejuvenating your skin,Strengthening your immune system, Enjoy the long list below...

Rejuvenates Skin

Improves Sleep

Recovers Faster

Burns Calories

Releases Toxins

Stimulates Happy Chemicals

Reduces Stress

Boosts Collagen

Reduces Puffiness

Eliminates free radicals

Promotes Blood Circulation

Boosts immunity

## Operating instruction

1. Connect the controller and connect the power.
2. Press the power button.
3. Press the photon button to turn on the photon function.
4. Press temperature ▲ ▼ adjust temperature.
5. Press time ▲ ▼ adjust the use time.
6. After use, please unplug the power.



## Care and storage:

1. Do not use sharp items with the product.
2. After Use, Clean the product and keep it avoiding direct sunlight.
3. Keep the product clean with a soft cloth.
4. Do not use water, Flammable or chemical solutions to clean the product. Do not store it in humid place.
5. Do not place heavy items on this product.
6. Do not immerse this product in water.